

31 July 1964  
GN:bb:366  
(997-112)

STATINTL

## TRIP REPORT

STATINTL Purpose: To Evaluate [REDACTED] Microdensitometers

STATINTL Persons

STATINTL Contacted: [REDACTED]

STATINTL  
STATINTL

STATINTL

STATINTL

During the day of 23 July, a series of tests<sup>1</sup> were conducted on a [REDACTED]  
Model III Microdensitometer. The instrument was operated by [REDACTED]

STATINTL

[REDACTED] indicated that [REDACTED] entered the field of microdensitom-  
etry in 1956 with the intention of making an instrument for in-house purposes.  
The Model III is a modified version of the original instrument.

The instrument has a scanning capability in only one direction. A Leitz  
stage is used to obtain sample rotation. The stage design is such that continuous  
scanning is limited to one inch in the x direction. Upon removal of the Leitz  
stage an eight inch scan can be accomplished, but rotation is not possible. The  
mensuration accuracy, as stated by [REDACTED] is  $\pm 15 \mu$ . The basic price  
for the Model III microdensitometer is [REDACTED]. A digitized system has recently  
been developed and is available at an additional cost.

STATINTL  
STATINTL  
STATINTL

A lower cost instrument, the [REDACTED] Model V, has been developed for  
commercial use. The Model V uses Quanta-log electronics, a Leitz stage and

STATINTL

Memorandum by [REDACTED] MM:bb:282, 22 June 1964.

Declass Review by NIMA/DOD

Approved For Release 2001/04/02 : CIA-RDP78B04747A000200010041-4

Trip Report  
Page 2

10 August 1964  
GN:bb:366 (997)

STATINTL Leitz optics. The instrument has the same scan limitation as the Model III. Its mensuration accuracy is estimated to be  $\pm 2-1/2 \mu$ . The overall cost is about [REDACTED]. We were informed that the price could be reduced considerably if a number of instruments were purchased.

There are no provisions on either the Model III or the Model V for handling roll film. A mirror system is used for viewing when the instrument is not in the scanning operation.

An automatic focusing device has been designed for use with their Model III microdensitometer. At the time of the visit the focusing device was under modification and would not be operational until the first part of August. Tentative arrangements were made for us to return when the system is operational. This will also give us an opportunity to discuss a number of items where information was initially lacking and which [REDACTED] is currently checking. STATINTL

STATINTL [REDACTED] indicated that [REDACTED] hopes to introduce a microdensitometer in the near future which would be designed specifically for commercial use. STATINTL He indicated the instrument would consist of a [REDACTED] Computer, STATINTL [REDACTED] electronics and [REDACTED] optics. The instrument would be assembled at [REDACTED] and special features would be available. STATINTL

[REDACTED]  
Ext. 547