

CONFIDENTIAL

October 2, 1964

[redacted]
Post Office Box 9642
Rosslyn Station
Arlington, Virginia 22209

25X1

Attention: [redacted]

25X1

Gentlemen:

[redacted] is pleased to submit for your evaluation a Cost and Technical Proposal for the design and manufacture of one (1) or two (2) pair of Low Power Wide Field Objectives for use with the [redacted] High Power Stereoviewer. Copies of our Technical Proposal are enclosed, as are copies of an outline drawing of the proposed objective.

25X1

25X1

A fixed price contract is anticipated. Our Cost Analysis is included as Attachment A.

Delivery of the end item is estimated at within 12 weeks after authorization to proceed.

If you have any questions concerning this proposal, please contact the writer directly.

25X1

GJJ:bp

[redacted]
Photogrammetric Contracts Section

CONFIDENTIAL

DDR-Dupe

Group 1
Excluded from automatic
downgrading and
declassification

Page Denied

25X1

TECHNICAL DESCRIPTION
and
DESIGN OBJECTIVES
for
PROTOTYPES
of a
[REDACTED]
WIDE FIELD OBJECTIVE

25X1

This lens is a low magnification, wide field microscope objective. It is approximately parfocal with the higher magnification objectives used with the [REDACTED] Dyna-Zoom Laboratory Microscope. This objective may be used with the four position nosepiece without interfering with the other objectives.

25X1

The real field of a microscope is equal to the field of the eyepiece divided by the magnification of the objective. The only way to get a significant increase in field is to lower the magnification of the objective. Lower magnification necessarily means increased working distance and below magnifications of approximately 3.5X the objectives cease to be parfocal.




The objective lens covered by this description has a prism cluster in the lower conjugate that folds up the optical path enabling the lens to have a shoulder-to-object distance the same as that of the other lenses used on the DynaZoom Microscope. This prism cluster does not invert or revert the image.

The nominal magnification of this lens is 1.3X. The field depends on the eyepiece. With the [REDACTED] 10X Wide Field Eyepiece, Catalog Number 31-05-60, the real field is approximately 14.0mm.

25X1

No resolution requirement is specified for the prototype lot since the primary objective is a wider field for initial location of areas of interest. Specifications for follow-on quantities will be determined after evaluation of the prototype.

25X1



ASS'Y SHOWN ON _____

C 588334

CHANGES

L
O
C
NUMBER
DATE

Z	
Y	
X	
W	
U	
T	
S	
R	
P	
N	
M	
L	
K	
H	
G	
F	
E	
D	
C	
B	
A	

DO NOT SCALE—WORK TO FIGURES
REPORT ERRORS—UNLESS OTHERWISE
SPECIFIED TOLERANCES ARE AS NOTED

TYPE FRAC. DEC. ANG.

MACHINED
CAST-MOLDED

THREADS—CL-2 UNIFIED

APPROVALS	SECT. HEAD
APPROVED FOR PROTOTYPE MANUFACTURING	
APPROVED FOR PRE-AUTHORIZATION TOOLING AND PARTS	
AUTHORIZED FOR PRODUCTION	

COSMETIC QUALITY - GRADE D UNLESS OTHERWISE SPECIFIED	MATL.	UNIT
	FINISH	SCALE 2 X SIZE
		DRWN
		CHKD
		APRD
		DATE 9-17-64

WIDE FIELD OBJ. 1.3 X

OUTLINE DRAWING

C 588334 -

Approved For Release 2006/12/11 : CIA-RDP78B04770A000100020007-3

Page Denied